

CENTRI[®]

A Fillauer Company

PEVA



TARSO® Shoe

Tarso® Shoe

The Tarso® Shoe is used for the Ponseti method of treating a club foot.

The shoe has an open toe with 2 threaded holes for connecting to the Dennis Browne bar.

- Now with deeper and rounder heel
- Padded and free standing tongue
- Leather lining
- Extra strengthening on the medial side



TARSO® Shoe		
Straight Medium Last		
Product no.	Size	Last Length
697 000 000	0000	74 mm
697 001 000	000	82 mm
697 001 100	00	90 mm
697 001 110	0	98 mm
697 001 111	1	107 mm
697 001 115	1,5	111 mm
697 001 112	2	117 mm
697 001 125	2,5	121 mm
697 001 113	3	125 mm
697 001 135	3,5	129 mm
697 001 114	4	134 mm
697 001 145	4,5	138 mm
697 001 150	5	142 mm
697 001 155	5,5	146 mm
697 001 116	6	151 mm
697 001 117	7	157 mm
697 001 118	8	166 mm
697 001 119	9	174 mm

TARSO® Shoes		
Straight Narrow Last		
Product no.	Size	Last Length
697 200 000	0000	74 mm
697 201 000	000	82 mm
697 201 100	00	90 mm
697 201 110	0	98 mm
697 201 111	1	107 mm
697 201 115	1,5	111 mm
697 201 112	2	117 mm
697 201 125	2,5	121 mm
697 201 113	3	125 mm
697 201 135	3,5	129 mm
697 201 114	4	134 mm
697 201 145	4,5	138 mm
697 201 150	5	142 mm
697 201 155	5,5	146 mm
697 201 116	6	151 mm
697 201 117	7	157 mm
697 201 118	8	166 mm

TARSO® Shoes		
Full Abduktion Medium Last		
Product no.	Size	Last Length
697 101 000	000	83 mm
697 101 100	00	91 mm
697 101 110	0	99 mm
697 101 111	1	108 mm
697 101 115	1,5	112 mm
697 101 112	2	118 mm
697 101 125	2,5	121 mm
697 101 113	3	125 mm
697 101 135	3,5	129 mm
697 101 114	4	133 mm
697 101 145	4,5	137 mm
697 101 150	5	141 mm

PEVA Shoes

PEVA Shoes are designed to be used to help maintain correction in babies with PEVA foot.

- Padded and free standing tongue
- Leather lining
- Extra strengthening of the heel cap



PEVA Shoes		
Straight Medium Last		
Product no.	Size	Last Length
697 800 012	12	73 mm
697 800 013	13	80 mm
697 800 014	14	86 mm
697 800 015	15	93 mm
697 800 016	16	100 mm
697 800 017	17	106 mm
697 800 018	18	113 mm
697 800 019	19	120 mm
697 800 020	20	126 mm
697 800 021	21	133 mm
697 800 022	22	140 mm
697 800 023	23	147 mm
697 800 024	24	153 mm

PEVA Shoes		
Straight Narrow Last		
Product no.	Size	Last Length
697 805 012	12	73 mm
697 805 013	13	80 mm
697 805 014	14	86 mm
697 805 015	15	93 mm
697 805 016	16	100 mm
697 805 017	17	106 mm
697 805 018	18	113 mm
697 805 019	19	120 mm
697 805 020	20	126 mm
697 805 021	21	133 mm
697 805 022	22	140 mm
697 805 023	23	147 mm
697 805 024	24	153 mm

Dennis Brown Bar

Dennis Browne Bar

Dennis Browne bar is an adjustable bar with foot plates.

The foot plates have 2 larger holes for the screws.

All plates have a scale with 10 degrees steps and are clearly marked.

Includes adjustable tools, allen wrench and screws.

- Adjustable bar
- Incl. adjustable tools, allen wrench and screws



Dennis Browne Bar		
Product no.	Length	Description
696 012 195	17,5 - 27,5 cm	Small plates (7,5 cm)
696 012 203	22,5 - 37,5 cm	Small plates (7,5 cm)
696 012 196	17,5 - 27,5 cm	Large plates (10 cm)
696 012 204	22,5 - 37,5 cm	Large plates (10 cm)

Horton Click Bar

Horton Click Bar

The Horton click is a foot rotation bar that allows the feet to be internally or externally rotated.

The Horton Click is easier to fit on the child because you fit the shoes on independently then click them on the foot rotation bar. This allows for better positioning of the foot in the shoe which reduces the child from being able to kick the shoes off. One of the many advantages of the Horton Click is the ease of dressing the child. With just a click you can quickly remove the bar from the shoes, change the child, and click the shoes back on the bar.



Horton Click Bar		
Product no.	Length	Description
697 000 011	20,3 - 30,5 cm	Quick release, Adjustable bar with 10 degree Kicker and Screws
697 000 012	15,2 - 20,3 cm	Quick release, Adjustable bar with 10 degree Kicker and Screws
697 000 013	20,3 - 30,5 cm	Quick release, Adjustable bar
697 000 014	15,2 - 20,3	Quick release, Adjustable bar



Dobbs Bar

Dobbs Bar

The Dobbs Bar lets children move and kick their legs independently, which gives them new found freedom when it comes to crawling and being more active. Children fuss less and sleep better in the Dobbs Bar



- Children can kick and move their legs independently so they can crawl and be more active.
- No more blisters and red spots on the back of the heels, because the brace moves with the child, friction on the heel is greatly reduced.
- The independent movement and greater mobility of the brace prevents the kids from escaping and pulling out of their boots.
- Children sleep better in the Dobbs Bar because they can get comfortable.
- Children can sleep with one leg straight and the other bent which reduces stress on the child's knees.
- The Dobbs Bar attaches to our Open-Toe Boots and is completely adjustable from 8 to 14 inches to allow for growth.

Dobbs Bars		
Product no.	Length	Description
697 000 010	20,3 - 30,4 cm	Plates (10 cm)

EZ CRAWL



- Allows greater freedom of movement for crawling and other activities
- Adheres to modern treatment protocol for Talipes Equinovarus (clubfoot)
- Quick and easy detachment for diaper changes, donning, and doffing
- Permits standing and full weight bearing
- Adjustable bar width to accommodate growth
- Allows footwear with orthosis
- Lightweight and versatile design
- Developed by leading clinic pediatric orthotists



EZ Crawl Bar	
Product no.	
	696 012 400



A Fillauer Company