

PROVIDENCE NIGHTTIME BRACE

Features & Benefits

Providence nighttime braces provide a viable alternative for many juvenile and adolescent patients. Recent studies indicate that this type of regimen can yield impressive initial correction as well as long-term success.

- Over corrective
- Nighttime use
- Controls opposite forces
- Controls directly lateral forces
- Controls rotational forces
- Produced through CAD/CAM
- Polypropylene material
- Special measurement board

Recent Study

In a recent study reported in *The Journal of Spine* (September 2001), 102 female patients with Risser 0, 1, and 2 were treated using only the Providence brace. Patients wore the brace between eight and ten hours a night for six or seven days out of the week, a far cry from the time requirements traditionally associated with full-time spinal bracing.

The brace achieved a 74% success rate.



PROVIDENCE NIGHTTIME BRACE

- Over corrective
- Nighttime use
- Controls opposite forces
- Controls both directly lateral and rotational forces
- Polypropylene material
- Special measurement board

Ordering Number	
Part no	Product description
688 000 101	Measurement Board
688 500 000	Providence Brace finished
688 500 001	Providence Brace unfinished

