

## Motionfoot



A revolutionary product providing unparalleled range of motion - 50° range of plantar/ dorsi flexion with a pending hydraulic ankle.

Prosthetist or wearer can adjust both directions of resistance independently, to accommodate a variety of tasks and activities.

The ankle can also be manually locked at neutral for "fixed ankle" activities, like driving, standing, etc.

- Independent hydraulic adjustments for plantar flexion resistance and dorsi flexion resistance
- Total ROM - 50°
- Ankle can be locked in neutral position
- Adjustable keel dynamics
- Foot sizes 24 - 30 cm
- Weight rating of 100 kg

## Ordering information

Product number	Description
732 001 220	Motion foot size 22 cm
732 001 230	Motion foot size 23 cm
732 001 240	Motion foot size 24 cm
732 001 250	Motion foot size 25 cm
732 001 260	Motion foot size 26 cm
732 001 270	Motion foot size 27 cm
732 001 280	Motion foot size 28 cm
732 001 290	Motion foot size 29 cm
732 001 300	Motion foot size 30 cm

## Motion foot



## Ordering information

Product number	Description
452 213 033	Foot shell, left size 22
452 313 033	Foot shell, left size 23
452 413 033	Foot shell, left size 24
452 513 033	Foot shell, left size 25
452 613 033	Foot shell, left size 26
452 713 033	Foot shell, left size 27
452 813 033	Foot shell, left size 28
452 913 033	Foot shell, left size 29
453 013 033	Foot shell, left size 30
462 213 033	Foot shell, right size 22
462 313 033	Foot shell, right size 23
462 413 033	Foot shell, right size 24
462 513 033	Foot shell, right size 25
462 613 033	Foot shell, right size 26
462 713 033	Foot shell, right size 27
462 813 033	Foot shell, right size 28
462 913 033	Foot shell, right size 29
463 013 033	Foot shell, right size 30

## Foot shell

